


HPLC Analysis of the Catecholamine Pathway



Primesep 200 separates catecholamines in the catecholamine pathway in 10 minutes. Phenylalanine, tyrosine, DOPA, dopamine, norepinephrine, and epinephrine are baseline resolved by a combination of reversed-phase, ion-exchange, and ion-exclusion mechanisms. Excellent peak shape results with a mass spec compatible mobile phase of water, acetonitrile (MeCN, ACN) and trifluoroacetic acid (TFA) with UV detection at 210 nm.

Method Parameters

Column	Primesep 200, 3.2x150 mm, 5 µm, 100 Å
Mobile Phase	MeCN/H ₂ O – 10/90%
Buffer	TFA – 0.1%
Flow Rate	0.5 mL/min
Detection	UV, 210 nm

Quelle: <https://sielc.com/Application-HPLC-Analysis-of-the-Catecholamine-Pathway>