




Effect of Buffer Concentration on Resolution of Neutral Compounds (Prednisolone and D-Glucose Pentaacetate)

<https://sielc.com/Application-Effect-of-Buffer-Concentration-on-Resolution-of-Neutral-Compounds-Prednisolone-and-D-Glucose-Pentaacetate>

Chromatogram



Column: Primesep P
Column size: 150 x 4.6 mm
Flow rate: 1.0 mL/min
Detection: ELSD



Description

Effect of buffer concentration on retention and selectivity of HPLC separation of prednisolone and glucose pentaacetate was studied on a mixed-mode Primesep P cation-exchange column.