

<https://sielc.com/hplc-method-for-analysis-of-nirmatrelvir-in-paxlovid-nirmatrelvir-tablets>

Column:	Primesep P
Column size:	3.2 × 100 mm, 5 µm
Column part number:	P-3.2.100.0510
Mobile phase:	Gradient MeCN – 10-90%, 10 min
Buffer:	No
Flow rate:	0.5 mL/min
Detection:	ESI Negative SIM

AugustaChrom – Dr. Roland Kießwetter · Langemarckstr. 21c · 86368 Gersthofen · Tel: 0821 599 783 70 · Fax: 0821 599 783 71 · info@hplc-shop.de · hplc-shop.de

Method Parameters

Mobile Phase	Gradient MeCN – 10-90%
Buffer	No
Flow Rate	0.5 ml/min
Detection	UV 210nm
Class of Compounds	Drug
Analyzing Compounds	Nirmatrelvir

HPLC Column Used

Primesep P, 3.2 x 100 mm, 5 µm, 100 A, dual ended

[Order this column at hplc-shop.de →](https://www.hplc-shop.de)